


3) Sample extraction and preparation (PRE-CEV INTERVIEW)

- Weight 500 mg of the chosen infuse directly in Falcon by using the analytical balance. Record the weight performed.
- Add to the sample 10 mL of MillQ water, stir for 1 minute, then let it extract statically (*Which kind of extraction is?*).
- Filter the extract by using a syringe filter (PTFE): (i) screw the filter to the syringe; (ii) place the syringe filter on the 50 mL falcon; (iii) load the extract in the syringe; (iv) perform the filtration.
- Perform the calculation to prepare the dilutions described below:

Dilute the extract 1:5, 1:10, 1:20, 1:50, 1:100, and 1:200 in a final volume of 1000 μL by using buffer pH 10 as solvent, using **serial dilution** methods.

See the figure and the suggestion below:

- Directly from the extract prepare serially the dilutions 1:5, 1:10, and 1:20
- From the solution 1:5, prepare the solution diluted 1:50
- From the solution 1:10, prepare the solution diluted 1:100
- From the solution 1:20, prepare the solution diluted 1:200



PAY ATTENTION: PREPARE THE SOLUTIONS AFTER THE CEV INTERVIEW.

Before starting: (i) understand the concentrations to be diluted serially; (ii) get the dilution factor among the concentration set and sign it in the figure; (iii) see the suggestion below; (iv) fill the table with the calculated volume considering the **final volume of 1 mL**.

DILUTIONS	1:5	1:10	1:20	1:50	1:100	1:200
	μL	μL	μL	μL	μL	μL
Extract						
Buffer pH 10						

4) External calibration curve building-up and sample analysis (POST-CEV INTERVIEW)

- Prepare the standard solutions according to **point 2)**. Mix each solution by pipetting to ensure a homogenous solution.
- Prepare the sample dilutions according to **point 3)**. Mix each solution by pipetting to ensure a homogenous solution.
- Align the $n= 17$ colorimetric strips on the white support including also the blank of reaction (*What does it mean?*).
- Onto each colorimetric strip load 100 μL of standard solutions (5 - 200 μM) and diluted extract (1:5-1:200) according to increasing concentration and dilution degree and let them react for 20 minutes at 50°C on the heated plate (up to strip drying). PAY ATTENTION! YOU SHOULD RUN THE BLANK OF REACTION (Buffer pH 10 only) TOGETHER TO THE STANDARDS AND SAMPLES.
- At the end of the reaction, recover the reacted colorimetric strips.
- Open the RGB color detector and collect pictures ensuring a good and uniform lighting condition.
- Using the RGB color detector record the Red (R) and Blue (B) hues for each reacted colorimetric strip. Before collecting the value adjust the pointer/sampler size (*Why we don't collect the Green hue?*).
- Report the collected data in an Excel file.

COMPARISON CHART FOR DIRECT QUANTIFICATION:

