

# Oocyte Quality Assessment by Brilliant Cresyl Blue (BCB) staining

## AIM OF THE STAINING:



Used to select more competent oocytes prior to IVM in various species, including pigs, mice, goats, cattle, and buffaloes

## HOW BCB WORKS:


The staining allows for the determination of the glucose-6-phosphate dehydrogenase (G6PDH) activity, which converts BCB dye from blue to colorless.

The oocytes, which are still in the growth phase, have a high G6PDH activity and show a colorless cytoplasm (BCB-). However, oocytes that have completed their growth have low levels of G6PDH and show a blue coloration of the cytoplasm (BCB+).

Research studies have demonstrated that BCB+ oocytes have a significantly higher blastocyst developmental rate than BCB- oocytes, suggesting that the quality of BCB+ oocytes is higher than that of BCB- oocytes.



BCB-



BCB+

## Oocyte Quality Assessment by Brilliant Cresyl Blue (BCB) staining

*BCB Molecular weight (MW): 385.96 g/mol*

- Prepare a BCB stock solution 1mg/mL, this correspond to a concentration of 2.6mM.

Let's verify:

$$g = MW \times V \times M \longrightarrow M = \frac{g}{MW \times V}$$

- To use the formula we need to express “mg as g” and “mL as L”:

$$1\text{mg/mL} = 0.001\text{g}/0.001\text{L}$$

- Let's put the values in the formula:

$$M = \frac{0.001 \text{ g}}{385.96 \times 0.001\text{L}} = 0.0026 \text{ M} \longrightarrow 2.6 \text{ mM}$$

- The BCB powder should be dissolved in 10mL of PBS at the concentration 1mg/mL to have the stock solution.

How many gr?

.....

- The medium for incubation of CoC with BCB is PBS + 0.5% BSA (mPBS). Prepare 10mL of this solution. How many gr of BSA?

.....

- BCB Working solution for pig CoC → 65 $\mu$ M

Calculate the dilution factor for BCB:

.....

How many ul of BCB stock solution per 10mL of mPBS (to have the BCB working solution)?

.....

### **WE ARE READY FOR STAINING CoC:**

1. Collect CoC in a multi-well plate
2. Wash twice CoC with mPBS (0.5mL each)
3. Move CoC into the BCB working solution and incubate 1h 38°C
4. Wash twice CoC with complete growth medium (alphaMEM + 10% FBS + 1%P/S + 1% Amphotericin)
5. Observe CoC with a stereomicroscope

# Oocyte nuclei staining with Propidium Iodide (PI)

## **AIM OF THE STAINING:**

Chromatin configuration characterization to distinguish fully grown oocytes from growing ones

## **PROTOCOL:**

1. Denude the oocytes.
2. Transfer them to PBS.
3. Fix in 1% PFA in PBS for 15 minutes.
4. Wash twice with PBS.
5. Permeabilize with 0.5% Triton X-100 in PBS for 20 minutes.
6. Wash twice with PBS.
7. Stain with 1 mg/mL Propidium Iodide for 30 minutes in the dark. The stock solution is 1 g/mL.
8. Wash twice with PBS.
9. Mount the stained oocytes on a slide and cover with a coverslip.
10. Observe using a fluorescence microscope or a confocal station.