



How to run Alpha assay:


- How to setup an Alpha assay
- Make your own assay!

➤ Samples:

- Phenol red and hemoglobin: choose AlphaLISA beads instead of AlphaScreen beads
- Pay attention to **RPMI** and **tissue homogenates**: they can contain biotin which will bind to streptavidin-coated beads
 - Biotin-free kits are available
 - Carefully plan order of addition (pre-incubate biotinylated antibody + streptavidin-beads)
- 1% FBS or 0.1% BSA can be added to cell culture medium to enhance assay sensitivity
- Serum samples: serum should not exceed 10% of final assay
- If necessary, **dilute your samples with AlphaLISA buffer**

➤ Standard curve:

- Prepare:
 - standard curve in the **same matrix** as samples (lysis buffer, serum without analyte) and
 - **standard curve in AlphaLISA buffer in parallel** to check interference from matrix.
- If samples are diluted in AlphaLISA Buffer, prepare a standard curve in **diluted matrix** (example, serum diluted in AlphaLISA buffer)




AlphaLISA beads contain europium:

- highly intense and spectrally defined signal
- high signal/background ratio
- works in different samples including serum, plasma & cell culture supernatant

Donor beads: effect of light exposure

700 lux light exposure



Lighting condition	From (lux)	To (lux)	Mean value (lux)
Pitch Black	0	10	5
Very Dark	11	50	30
Dark Indoors	51	200	125
Dim Indoors	201	400	300
Normal Indoors	401	1000	700
Bright Indoors	1001	5000	3000
Dim Outdoors	5001	10,000	7500
Cloudy Outdoors	10,001	30,000	20,000
Direct Sunlight	30,001	100,000	65,000



- Switch off lights!
- Avoid direct day light (≤ 100 lux)
- Protect Donor beads vials with aluminum
- Cover plate (with another plate)
- Incubate plate in the dark (drawer)
- Use green filters



Alpha
 Absorbance
 Fluorescence top & bottom
 Ultra-sensitive luminescence
 TRF
 Label-free

Please select a protocol and press the 'Run' button above

Valid	Technology	Name
	Abs	Protein Concentration Absorbance
	Abs	Spectrum Scan Absorbance
	Abs	Nucleotide Purity Absorbance
	Alpha	Alpha 96-well
	Alpha	Alpha 384-well Low Volume
	Alpha	Alpha 384-well
	Alpha	Performance IPA Alpha

	1	2	3	4	5	6	7	8	9	10	11	12
A	STD1	STD2	STD3	STD4	STD5	STD6	STD7	STD8	STD9	STD10	BL	BL
B	STD1	STD2	STD3	STD4	STD5	STD6	STD7	STD8	STD9	STD10	BL	BL

Curve fitting UNK calculation

Source: Meas A (Alpha)

Fitting: 4PL

Weighted: Weighted

X-axis scale: Logarithmic

Y-axis scale: Logarithmic


Concentration unit: pg/ml

Advanced options


STD	Concentration
1	30000
2	10000
3	3000
4	1000
5	300
6	100
7	30
8	10
9	3
10	1

- ▶ Plot signal vs concentration of analyte (Log scale for either or both axes)
- ▶ Analyze data according to a nonlinear regression using the **4-Parameter Logistic** equation (Sigmoidal Dose-Response curve with variable slope)
- ▶ Add $1/Y^2$ weighting
- ▶ Values at maximal concentrations of analyte after the hook point must be **removed**
- ▶ **LDL (12 blk) = average + 3 x st.dev.**

St. curved in Alpha buffer



St. curve in diluted serum




Notice the matrix effect of serum




Make your own assay


Develop your own assay with Alpha Toolbox

- ▶ Beads can be coated with antibodies or other binding molecules to develop virtually any immunoassay




- ▶ Virtually any assay can be developed, as long as you can bring the beads together





 primary Ab



 Protein A (or secondary IgG)


Direct conjugation

- Simple reaction (1 night incubation)
- Good control on binding events
- Random coupling (antibodies aligned randomly)
- Simple stoichiometry in Alpha assay

Indirect binding

- No custom bead conjugation required
- Save expensive primary Ab
- Antibodies aligned optimally
- **Additional equilibrium introduced**
- **Needs high affinity primary Ab + beads with selective secondary Ab or proteins A/G/L**

Reductive amination:



- ▶ Monoclonal or purified polyclonal antibodies will perform best; **avoid anti-sera**
- ▶ Antibody concentration must be at least **1 mg/ml**
- ▶ **Antibody must not be in any amine-based buffer**, including TRIS, glycine, bicine, tricine, etc... Perform a buffer exchange in PBS pH 7.4
- ▶ Antibody solution **must not contain any protein (such BSA or gelatin) and glycerol**

Antibody assesment and pre-treatment (if needed)

Wash raw beads

20 min.

Assemble conjugation reaction in a vial

10 min.

Overnight incubation

(18-24 h 37°C)

Carboxymethoxylamine (CMO) blocking step

60 min. 37°C

Washing steps

60 min.

Many types of derivatized beads

Donor beads:

Streptavidin
Nickel chelate (His-tag)
GSH
Protein A
Anti-FLAG
Anti-DIG
Anti-mouse IgG
Anti-Rabbit IgG
Strep-Tactin
LCA (Lens Culinaris Agglutinin)

Acceptor beads:

Antibody capture

Protein A
Protein G
Protein L
Anti-human IgG1
Anti-human IgG4
Anti-human IgG
Anti-rabbit IgG
Anti-mouse IgG
Anti-rat IgG
Anti-goat IgG
Anti-Mouse IgM
Anti-Chicken IgY
Anti-Sheep IgG
Anti-bovine IgGA, IgG, IgG1, IgG2, IgGM

Anti-mouse:
IgE
IgM
IgG1
IgG2a
IgG2b
IgG3
(isotyping)

Fusion Tag detection

Nickel chelate (His-tags)
Anti-His
Glutathione (GSH)
Anti-GST
Anti-c-myc
Anti-FLAG
Anti-DIG
Anti- HA
Anti-FITC
Anti-V5
Anti-GFP
Anti-MBP
LCA
Strep-Tactin
Streptavidin

Streptavidin donor beads are generally preferred



- ▶ Protein A interacts strongly with:
 - Human IgG1, IgG2 and IgG4 and total IgG
 - Mouse IgG2A and IgG2B and total IgG
 - Rabbit total IgG

- ▶ Protein G binds to all subclasses of human and mouse IgG and to rat, goat, sheep, guinea pig, rabbit, cow, pig and horse antibodies


- ▶ Protein L binds to:
 - Total human IgG, IgM, IgA, IgE, IgD
 - Mouse IgG
 - Rat IgG
 - Binds poorly to mouse IgM and rabbit IgM
 - Does not bind to rabbit, sheep, goat and bovine IgG and IgM

Titration of biotinylated antibody




The conc. of biotinylated antibody must be titrated to avoid Hook effect




Hook effect in homogeneous assays




The conc. of biotinylated antibody must be titrated to avoid Hook effect



Choose antibody conc. giving strong Alpha signal, without Hook effect



 primary Ab



 Protein A (or secondary IgG)

Direct conjugation

- Simple reaction (1 night incubation)
- Good control on binding events
- Random coupling (antibodies aligned randomly)
- Simple stoichiometry in Alpha assay

Indirect binding

- No custom bead conjugation required
- Save expensive primary Ab
- Antibodies aligned optimally
- **Additional equilibrium introduced**
- **Needs high affinity primary Ab + beads with selective secondary Ab or proteins A/G/L**

Check beads compatibility (antibody cross-reactivity):

ACCEPTOR BEADS	DONOR BEADS												
	Streptavidin (6750002)	Anti-FLAG (AS103)	Anti-Mouse (AS104)	Anti-Rabbit (AS106)	Streptactin (AS108)	Ni chelate (AS101)	Glutathione (6766300)	Protein A (AS102)	Anti-DNP (AS111)	Anti-HRP (AS109)	Anti-Rat (AS110)	Amb-DIG (AS108)	Anti-GFP (AS112)
Streptavidin (AL125)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Protein L (AL126)	✓	✓	✓	✓	✓	✗	✓	✓	✓	✓	✗	✓	✗
Anti-FITC (AL127)	✓	✓	✗	✓	✓	✓	✓	✗	✓	✓	✗	✗	✓
Anti-His (AL128)	✓	✓	✗	✓	✓	✓	✓	✗	✓	✓	✗	✗	✓
Anti-V5 (AL129)	✓	✓	✗	✓	✓	✓	✓	✗	✓	✓	✗	✗	✓
Anti-Mouse IgM (AL130)	✓	✓	✗	✓	✓	✓	✓	✗	✓	✓	✗	✗	✓
Anti-Chicken IgY (AL131)	✓	✓	✓	✗	✓	✓	✓	✗	✓	✓	✗	✗	✓
Anti-Sheep IgG (AL132)	✓	✓	✓	✗	✓	✓	✓	✗	✗	✓	✗	✗	✓
Anti-GFP (AL133)	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✗	✓	✓
Anti-MBP (AL134)	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✗	✓	✓
Streptactin (AL136)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✓
Protein A (AL101)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✓
Protein G (AL102)	✓	✓	✗	✗	✓	✓	✓	✗	✓	✓	✗	✗	✓
Anti-Human IgG (AL103)	✓	✓	✗	✓	✓	✓	✓	✗	✓	✓	✗	✗	✓
Anti-Rabbit IgG (AL104)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✓
Anti-Mouse IgG (AL105)	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✓
Anti-Rat IgG (AL106)	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✓
Anti-Goat IgG (AL107)	✓	✓	✓	✗	✓	✓	✓	✗	✓	✓	✗	✗	✓
Ni chelate (AL108)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Glutathione (AL109)	✓	✓	✓	✓	✓	✓	✓	✗	✓	✓	✗	✗	✓
Anti-GST (AL110)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Anti-cMyc (AL111)	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✗	✓	✓
Anti-FLAG (AL112)	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✗	✓	✓
Anti-Digoxigenin (AL113)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✗	✓	✓
Anti-DNP (AL173)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Anti-HRP (AL171)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Anti-Protein C tag (AL172)	✓	✓	✗	✓	✓	✓	✓	✗	✓	✓	✗	✗	✓
Anti-HA (AL170)	✓	✓	✓	✗	✓	✓	✓	✗	✓	✓	✗	✗	✓
Anti-Mouse IgG (AL164)	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✓
Anti-Human IgG1 (AL153)	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✓
Anti-Human IgG2 (AL154)	✓	✓	✗	✗	✓	✓	✓	✓	✓	✓	✗	✗	✓
Anti-Human IgG4 (AL156)	✓	✗	✗	✓	✓	✓	✓	✓	✓	✓	✗	✗	✓
Anti-Mouse IgE (AL161)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Anti-Mouse IgG1 (AL157)	✓	✗	✗	✓	✓	✓	✓	✓	✓	✓	✗	✗	✓
Anti-Mouse IgG2a (AL158)	✓	✗	✗	✓	✓	✓	✓	✓	✓	✓	✗	✗	✓
Anti-Mouse IgG2b (AL159)	✓	✓	✗	✓	✓	✓	✓	✓	✓	✓	✗	✗	✓
Anti-Mouse IgG3 (AL160)	✓	✓	✓	✗	✓	✓	✓	✗	✓	✓	✓	✓	✓
Anti-Mouse IgM (AL162)	✓	✗	✓	✓	✓	✓	✓	✓	✓	✓	✗	✗	✓
LCA (AL140)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓






Plate map for antibodies cross-titration

Choose antibodies conc. giving strong Alpha signal, without Hook effect



- ▶ Bead concentration: 10 - 40 $\mu\text{g/ml}$, keeping the ratio antibody/bead constant
- ▶ Incubation time: 30' – overnight
- ▶ Assay volume and sample volume
- ▶ Assay buffers
 - Tris, HEPES, PBS at different pH values
 - Mild detergents like Tween-20, CHAPS and Triton X-100 (0.01% - 1%)
 - Protein blockers, such as casein or BSA, at a concentration between 0.01 to 1%
 - Dextran 500: for some serum or plasma samples, it is important to add Dextran 500 at 1 mg/mL in order to prevent non-specific bead aggregation

More optimizations (2): order of addition

